

reconstruction of rfb files, Theo van Soest/UMC Utrecht, t.l.vansoest@umcutrecht.nl
based on RFX2PTW/Bert van der Leije/DDHK Rotterdam and Radpy/Stephen Flounder

++++++

syntax information:

*variables for structural information

UPPERCASE_WORD: variable with value to set conditions

*assignment

UPPERCASE_WORD:= value

*conditional blocks:

---if (condition)

..

*repeated structure:

<<MEASUREMENT REPEAT LOOP START>>

..
N:= value

<<MEASUREMENT REPEAT LOOP UNTIL N=0>>

*comments on structure

==DATA STRUCTURE START==

*used data types

smallint	signed integer 2 bytes (C++: int16_t)
double	double precision float 8 bytes
string	length-prefixed strings (Pascal-type string, size stored in byte[0])
string2	length-prefixed string with data starting at byte[2], byte[0]=length, byte[1]=0
[]	byte at current position of file pointer

*example: the structure starts with a string which takes 15 bytes,
followed by the number of groups as two-byte smallint

*type	length	contents	description
string	1+14	'Version:6.3.14'	VERSION:= 63; version string
smallint	2	70	GROUPS:= 70

++++++

==DATA STRUCTURE START==

*type	length	contents	description
READ_SUBHEADER:= true			
string	1+14	'Version:6.3.14'	VERSION:= 63; version string
smallint	2	70	GROUPS:= 70
smallint	2	-1	
string2	1+1+5	\$0'CBeam'	[number of bytes x0 CBeam] look for this structure

<<GROUP REPEAT LOOP START>>

```

string      1+19    'U02'          linac
?           2       ?              ?
double      8       10             energy [MeV]
smallint    2       0              modality (Photon=0,Electron=1,Proton=2,Neutron=3,Cobalt=4,Isotope=5)

---if (VERSION < v60)
?           2       ?              ?
double      8       ?              TG fieldsize [mm]
?           2       ?              ?
double      8       ?              AB fieldsize [mm]
---

smallint    2       -1              $3c
?           2       1              wedge type (0:Hard,1:Dynamic,2:Enhanced,3:Virtual,4:Soft)
smallint    2       0              ?
?           2       1              wedge angle
smallint    2       0              ?
?           2       1              gantry angle
smallint    2       90             ?
?           2       1              colimator angle
double      8       1000           ?
?           2       1              ssd [mm]
double      8       1000           ?
string      1+9     'Undefined'    sad [mm]
                                          APPLICATOR

---if (VERSION > v50) and ((APPLICATOR='') or (APPLICATOR[1]<>#0))
smallint    2       0              medium (0:air,1:water,2:film)
---

---if (VERSION > v53)
string      1+13    'Radiotherapie' clinic                                $6a
string      1+3     'UMC'          adres
string      1+0     ''             telefoon
string      1+0     ''             e-mail
?           2       1              ?
double      8       -200           TGmin field setting [mm]
?           2       1              ?
double      8       200            TGplus field setting [mm]
?           2       1              ?
double      8       -200           ABmin field setting [mm]
?           2       1              ?
double      8       200            ABplus field setting [mm]
smallint    2       0              gantryscale (0:CW_180_Down,1:CCW_180_Down,2:CW_180_Up,3:CCW_180_Up)
---

READ_SUBHEADER:= TRUE

<<MEASUREMENT REPEAT LOOP START>>

==SUBHEADER START: USE IF READ_SUBHEADER=TRUE==                                $a8

---if (READ_SUBHEADER = true)
---while ([ ]<>$80) and (([ ]<>$FF) or ([+1]<>$FF))

```

```

byte          1
---

--- ([]=$80)
?             1          $80          padding
---

--- ([]<=$80)
?             4          ?            ?
string2       1+1+13    x0'CProfileCurve' MEASUREMENTTYPE
---
---

==SUBHEADER END==

===measurement start==
longword      4          (04/02/2002 14:19) measurement datetime 7-8 hour difference for Europe, where
defined???
longword      4          (04/02/2002 14:19) modification datetime
byte          1          2          quantity (1:Relative optical density,2:Relative dose,3:Relative
ionisation,4:Absolute dose,5:Charge)

---if (VERSION > v50)
double        8          0          radius *0.1 [mm]
double        8          0          calibration factor

---if (VERSION > v51)
double        8          0          temperature
double        8          0          pressure
string        1+0        ''        calibration date
double        8          -1        detector offset [mm]
---

---if (VERSION <= v50)
double        8          ?          ?
double        8          ?          ?
---

string        1+4        'IC04'        detector name
smallint      2          4          detector type (1:Single diode,2:LDA-11,3:LDA-25,4:Ion chamber
(cylindrical),5:Ion chamber (plane parallel),6:Stereotactic,7:Film,8:CA24,9:BIS-2G')
string        1+7        'RM, Jdk'        operator
string        1+0        ''        comment

---if (quantity=4) or (measurement type='CSinglePointCurve') {absolute dose measurement | untested!!!}
?             6          ?          ?
smallint      2          n          number of measurements (NM)

```

?	2	?	?
---if (VERSION > v60)			
?	291	?	?

double	8	value	crossline
double	8	value	inline
double	8	value	beam
double	8	value	crossline_ref
double	8	value	inline_ref
---if (VERSION > v60)			
double	8	value	crossline_ref

double	8	value	beam_ref
---if (VERSION > v60)			
double	8	value	inline_ref

double	8	0	temperature
double	8	0	pressure
double * NM	8 * NM	values	all measurments
--- ([] <> #0) {zero value is error}			
?	2	?	

--- (end absolute)			
---if (VERSION > v53)			
smallint	2	1	mapping crossline: -z=-3, -y=-2, -x=-1, +x=1, +y=2, +z=3
smallint	2	-2	mapping inline : -z=-3, -y=-2, -x=-1, +x=1, +y=2, +z=3
smallint	2	-3	mapping beam : -z=-3, -y=-2, -x=-1, +x=1, +y=2, +z=3
smallint	2	1	measurements per point
smallint	1	1	?
double	8	10.5	scanspeed [mm/s]
smallint	2	13	?
smallint	2	6	?
smallint	2	2438	origin x (*0.1 [mm]) \$116
smallint	2	2425	origin y all at bottom corner of tank
smallint	2	3620	origin z center of tank is (2400,2400,2400)
double	8	0	isocenter crossline [mm]
double	8	0	isocenter inline [mm]
double	8	-1	isocenter beam [mm]
double	8	0	normalization crossline [mm]
double	8	2.5	normalization inline [mm]
double	8	121.4	normalization beam [mm]
double	8	100	norm value field

double	8	0	norm value reference	\$154
double	8	0	dark current field	\$15c
double	8	-96.3	dark current reference	
double	8	300	HV field	
double	8	300	HV reference	
smallint	2	42	gain field	\$17c
smallint	2	45	gain reference	\$17e
string	1+4	'HIGH'	range field	
string	1+4	'HIGH'	range reference	
double	8	94	water surface correction (motor unit) [mm]	
smallint	2	2	?	
smallint	2	1	?	
double	8	0	reference min	
double	8	0	reference max	
double	8	0	reference avg	
double	8	-1	?	
smallint	2	4	?	
double	8	1.09051254089422	renormalization value	
double	8	0	curve offset, equal to detector offset	
string	1+14	#13#10'commissioning'	setup comment	
smallint	2	0	?	
double	8	-243.4	point a: crossline [mm]	\$1d9
double	8	242.4	point a: inline [mm]	
double	8	18.9	point a: beam [mm]	
double	8	-243.5	point b: crossline [mm]	
double	8	-237.5	point b: inline [mm]	
double	8	19.4	point b: beam [mm]	
double	8	236.3	point c: crossline [mm]	
double	8	-237.5	point c: inline [mm]	
double	8	18.9	point c: beam [mm]	
double	8	236.4	point d: crossline [mm]	
double	8	242.4	point d: inline [mm]	
double	8	18.9	point d: beam [mm]	
?	10	?	?	

---if (VERSION <= v53)

---if (VERSION < v52)
 ? 8 ?

---if (VERSION >= v52)
 ? 10 ?

double	8	236.4	crossline start [0.1 mm]	\$243
double	8	0	inline start [0.1 mm]	
double	8	50	beam start [0.1 mm]	

double	8	-243.2	crossline end	[0.1 mm]	\$25B
double	8	0	inline end	[0.1 mm]	
double	8	50	beam end	[0.1 mm]	
smallint	2	1200	number of points (NP)		
<<REPEAT START NP times>>					
double	8	-243.2 (first point)	position	[0.1 mm]	
double	8	5.8 (first point)	value		
<<REPEAT END NP times>>					
byte	1	var	N		\$4d75

```

---if (N>0)
---while ([]<>$80) and (([]<>$FF) or ([+1]<>$FF))
  byte 1
---
```

```

---if ([]=$FF)
  READ_SUBHEADER:= TRUE
---
```

```

---if ([]<>$FF)
  READ_SUBHEADER:= FALSE
  byte 1
---
```

```

---if (N=0)
  READ_SUBHEADER:= TRUE
```

---if (VERSION>50)				
byte * 3	3	*	?	
byte	1	3	N	\$4d75

```

---if (N>0) and (VERSION>63) and (MEASUREMENTTYPE='CTMRCurve')
  READ_SUBHEADER:= FALSE
---
```

```

NEXT_MEASUREMENT
---
```

```
<<MEASUREMENT REPEAT LOOP UNTIL N=0>>
```

```
GROUPS:= GROUPS-1
```

```

---if (GROUPS>0)
--- WHILE ([]<>$80)
  byte 1
---
```

```
byte      1
---

<<GROUP REPEAT LOOP UNTIL GROUPS=0>>

<<GROUP REPEAT LOOP UNTIL GROUPS=0>>
==DATA STRUCTURE END==
```